


MDWF 1050

STERILIZING INSTRUMENTS

Presentation by Stephanie Duran



Step 1: Watch your hands.





Step 2: Put on Gloves and put on a mask

Step 3: Add the instrument in a bowl and the Cleaner.





Step 4: Scrub the instrument

Step 5: Rinse the Instrument





Step 6: Place the instrument into 10% solution chlorine, let it sit for 20 seconds.

Step 7: Place the instrument into another clean bowl and add lubricate spray.





Step 8: Boil distilled water into a clean deep pot. When the water is boiling, add the instrument for 20 minutes.

Step 9: After 20 minutes, remove the instrument with sterilized tongs.



Step 10: Place the instrument into a clean surface. Use the instrument immediately.

